

# PCSK9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01071-3**

## Capture Antibody Information

**Catalog Number:**  
84172-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241395G1  
**Reactivity:**  
human  
**GenBank:**  
NM\_174936.4

**Conjugate:**  
Unconjugated  
**Full name:**  
proprotein convertase subtilisin/kexin type 9  
**Gene ID:**  
255738

## Detection Antibody Information

**Catalog Number:**  
84172-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241395E11  
**Reactivity:**  
human  
**GenBank:**  
NM\_174936.4

**Conjugate:**  
Unconjugated  
**Full name:**  
proprotein convertase subtilisin/kexin type 9  
**Gene ID:**  
255738

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01071-3 targets PCSK9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCSK9 Recombinant antibody, PBS Only (Capture) 84172-4-PBS (241395G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCSK9 Recombinant antibody, PBS Only (Detector) 84172-3-PBS (241395E11). 100 µg. Concentration 1 mg/mL.

Alternative PCSK9 matched antibody pairs: MP00059-1, MP00059-2, MP00059-3, MP01071-1, MP01071-2, MP01071-4, MP50055-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01071-3, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-4-PBS. Detection antibody: 84172-3-PBS. Standard: Eg0362. Range: 1.563-100 ng/mL.