

# Kininogen 1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01081-1**

## Capture Antibody Information

Catalog Number: 84175-4-PBS	Clone ID: 241397G9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: kininogen 1
Isotype: IgG	GenBank: NM_001102416.3	Gene ID: 3827
Purification Method: Protein A purification	Immunogen Catalog Number: Eg1162	

## Detection Antibody Information

Catalog Number: 84175-3-PBS	Clone ID: 241397C8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: kininogen 1
Isotype: IgG	GenBank: NM_001102416.3	Gene ID: 3827
Purification Method: Protein A purification	Immunogen Catalog Number: Eg1162	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP01081-1 targets Kininogen 1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Kininogen 1 Recombinant antibody, PBS Only (Capture) 84175-4-PBS (241397G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: Kininogen 1 Recombinant antibody, PBS Only 84175-3-PBS (241397C8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01081-1, Kininogen-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84175-4-PBS. Detection antibody: 84175-3-PBS. Standard: Eg1162. Range: 0.781-100 ng/mL.