

For Research Use Only

# Mouse CD6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01099-3

## Capture Antibody Information

Catalog Number:  
84166-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241471A7

Reactivity:  
mouse

GenBank:  
NM\_009852.3

Immunogen Catalog Number:  
Eg1401

Conjugate:  
Unconjugated

Full name:  
CD6 antigen

Gene ID:  
12511

## Detection Antibody Information

Catalog Number:  
84166-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241471D4

Reactivity:  
mouse

GenBank:  
NM\_009852.3

Immunogen Catalog Number:  
Eg1401

Conjugate:  
Unconjugated

Full name:  
CD6 antigen

Gene ID:  
12511

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01099-3 targets CD6 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Mouse CD6 Recombinant antibody, PBS Only (Capture) 84166-4-PBS (241471A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD6 Recombinant antibody, PBS Only (Capture/Detector) 84166-2-PBS (241471D4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01099-3, Mouse CD6 Recombinant Matched Antibody Pair - PBS only. 84166-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1401. 84166-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



Serum of twelve mice was measured. The mouse CD6 concentration of detected samples was determined to be 2.83 ng/mL with a range of 1.10 - 4.98 ng/mL.