

COMMD9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01110-2

Capture Antibody Information

Catalog Number: 84251-2-PBS	Clone ID: 241505B11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: COMM domain containing 9
Isotype: IgG	GenBank: BC010892	Gene ID: 29099
Purification Method: Protein A purification	Immunogen Catalog Number: Ag35078	

Detection Antibody Information

Catalog Number: 84251-1-PBS	Clone ID: 241505A10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: COMM domain containing 9
Isotype: IgG	GenBank: BC010892	Gene ID: 29099
Purification Method: Protein A purification	Immunogen Catalog Number: Ag35078	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01110-2 targets COMMD9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: COMMD9 Recombinant antibody, PBS Only (Capture) 84251-2-PBS (241505B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: COMMD9 Recombinant antibody, PBS Only (Detector) 84251-1-PBS (241505A10). 100 µg. Concentration 1 mg/mL.

Alternative COMMD9 matched antibody pairs: MP01110-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01110-2, COMMD9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84251-2-PBS. Detection antibody: 84251-1-PBS. Standard: Ag35078. Range: 0.625-40 ng/mL.