

GCGR Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01118-2**

Capture Antibody Information

Catalog Number: 84237-1-PBS	Clone ID: 241501A6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: glucagon receptor
Isotype: IgG	GenBank: BC104854	Gene ID: 2642
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24861	

Detection Antibody Information

Catalog Number: 84237-4-PBS	Clone ID: 241501H8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: glucagon receptor
Isotype: IgG	GenBank: BC104854	Gene ID: 2642
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24861	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01118-2 targets GCGR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GCGR Recombinant antibody, PBS Only (Capture) 84237-1-PBS (241501A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: GCGR Recombinant antibody, PBS Only (Detector) 84237-4-PBS (241501H8). 100 µg. Concentration 1 mg/mL.

Alternative GCGR matched antibody pairs: MP01118-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01118-2, GCGR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84237-1-PBS. Detection antibody: 84237-4-PBS. Standard: Ag24861. Range: 0.313-20 ng/mL