

MLH1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01140-2**

Capture Antibody Information

Catalog Number: 84208-4-PBS	Clone ID: 241472H5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli)
Isotype: IgG	GenBank: BC006850	Gene ID: 4292
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2319	

Detection Antibody Information

Catalog Number: 84208-3-PBS	Clone ID: 241472D7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli)
Isotype: IgG	GenBank: BC006850	Gene ID: 4292
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2319	

Applications

Tested Applications: Cytometric bead array, Sandwich ELISA	Range: 0.313-40 ng/mL (Cytometric Bead Array) 62.5-4000 pg/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01140-2 targets MLH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: MLH1 Recombinant antibody, PBS Only (Capture) 84208-4-PBS (241472H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: MLH1 Recombinant antibody, PBS Only (Capture/Detector) 84208-3-PBS (241472D7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only


Selected Validation Data



Sandwich ELISA standard curve of MP01140-2, Human MLH1 Recombinant Matched Antibody Pair - PBS only. 84208-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2319. 84208-3-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL



The mean MLH1 concentration was determined to be 9.42 ng/mL in Jurkat cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP01140-2, MLH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84208-4-PBS. Detection antibody: 84208-3-PBS. Standard: Ag2319. Range: 0.313-40 ng/mL