

S100B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01165-1**

Capture Antibody Information

Catalog Number: 82271-12-PBS	Clone ID: 230169D2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: S100 calcium binding protein B
Isotype: IgG	GenBank: BC001766	Gene ID: 6285
Purification Method: Protein A purification	Immunogen Catalog Number: Ag7440	

Detection Antibody Information

Catalog Number: 82271-11-PBS	Clone ID: 230169B2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: S100 calcium binding protein B
Isotype: IgG	GenBank: BC001766	Gene ID: 6285
Purification Method: Protein A purification	Immunogen Catalog Number: Ag7440	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01165-1 targets S100B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: S100B Recombinant antibody, PBS Only (Capture) 82271-12-PBS (230169D2). 100 µg. Concentration 1 mg/mL.

Detection antibody: S100B Recombinant antibody, PBS Only (Detector) 82271-11-PBS (230169B2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01165-1, S100B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82271-12-PBS. Detection antibody: 82271-11-PBS. Standard: Ag7440. Range: 1.25-20 ng/mL.