

# ADAM17 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01169-1

## Capture Antibody Information

**Catalog Number:**  
84292-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241305F10  
**Reactivity:**  
human  
**GenBank:**  
BC136783  
**Immunogen Catalog Number:**  
Ag32418

**Conjugate:**  
Unconjugated  
**Full name:**  
ADAM metallopeptidase domain 17  
**Gene ID:**  
6868

## Detection Antibody Information

**Catalog Number:**  
84292-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241305D1  
**Reactivity:**  
human  
**GenBank:**  
BC136783  
**Immunogen Catalog Number:**  
Ag32418

**Conjugate:**  
Unconjugated  
**Full name:**  
ADAM metallopeptidase domain 17  
**Gene ID:**  
6868

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01169-1 targets ADAM17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAM17 Recombinant antibody, PBS Only (Capture) 84292-2-PBS (241305F10). 100 µg. Concentration 1 mg/mL.

Detection antibody: ADAM17 Recombinant antibody, PBS Only (Detector) 84292-1-PBS (241305D1). 100 µg. Concentration 1 mg/mL.

Alternative ADAM17 matched antibody pairs: MP01169-2, MP50392-1, MP50392-2, MP50392-3, MP50392-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01169-1, ADAM17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84292-2-PBS. Detection antibody: 84292-1-PBS. Standard: Ag32418. Range: 0.313-40 ng/mL.