

C7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01179-1**

Capture Antibody Information

Catalog Number: 84263-1-PBS	Clone ID: 241551E11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: complement component 7
Isotype: IgG	GenBank: NM_000587.4	Gene ID: 730
Purification Method: Protein A purification	Immunogen Catalog Number: Eg1202	

Detection Antibody Information

Catalog Number: 84263-5-PBS	Clone ID: 241551D6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: complement component 7
Isotype: IgG	GenBank: NM_000587.4	Gene ID: 730
Purification Method: Protein A purification	Immunogen Catalog Number: Eg1202	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01179-1 targets C7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Complement component C7/C7 Recombinant antibody, PBS Only 84263-1-PBS (241551E11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: C7 Recombinant antibody, PBS Only 84263-5-PBS (241551D6). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01179-1, C7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84263-1-PBS. Detection antibody: 84263-5-PBS. Standard: Eg1202. Range: 0.625-80 ng/mL