

RINT1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01190-3**

Capture Antibody Information

| | | |
|--|-------------------------------------|----------------------------------|
| Catalog Number: 84290-1-PBS | Clone ID: 241592E4 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: RAD50 interactor 1 |
| Isotype: IgG | GenBank: BC068483 | Gene ID: 60561 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag6145 | |

Detection Antibody Information

| | | |
|--|-------------------------------------|----------------------------------|
| Catalog Number: 84290-2-PBS | Clone ID: 241592F8 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: RAD50 interactor 1 |
| Isotype: IgG | GenBank: BC068483 | Gene ID: 60561 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag6145 | |

Applications

| | | |
|---|--|--|
| Tested Applications: Cytometric bead array | Range: 0.625-80 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|---|--|--|

Product Information

MP01190-3 targets RINT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RINT1 Recombinant antibody, PBS Only (Capture) 84290-1-PBS (241592E4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RINT1 Recombinant antibody, PBS Only (Detector) 84290-2-PBS (241592F8). 100 μ g. Concentration 1 mg/mL.

Alternative RINT1 matched antibody pairs: MP01190-1, MP01190-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01190-3, RINT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84290-1-PBS. Detection antibody: 84290-2-PBS. Standard: Ag6145. Range: 0.625-80 ng/mL