

PFKM Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01204-2**

Capture Antibody Information

Catalog Number:
84281-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241598F7

Reactivity:
human

GenBank:
BC021203

Immunogen Catalog Number:
Ag30840

Conjugate:
Unconjugated

Full name:
phosphofructokinase, muscle

Gene ID:
5213

Detection Antibody Information

Catalog Number:
84281-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241598A7

Reactivity:
human

GenBank:
BC021203

Immunogen Catalog Number:
Ag30840

Conjugate:
Unconjugated

Full name:
phosphofructokinase, muscle

Gene ID:
5213

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01204-2 targets PFKM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PFKM Recombinant antibody, PBS Only (Capture) 84281-1-PBS (241598F7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PFKM Recombinant antibody, PBS Only (Capture/Detector) 84281-2-PBS (241598A7). 100 µg. Concentration 1 mg/mL.

Alternative PFKM matched antibody pairs: MP01204-1, MP01204-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01204-2, PFKM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84281-1-PBS. Detection antibody: 84281-2-PBS. Standard: Ag30840. Range: 0.156-20 ng/mL