

Rat IL-1 beta Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01210-3**

Capture Antibody Information

Catalog Number:
84274-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241001G6
Reactivity:
rat
GenBank:
NM_031512

Conjugate:
Unconjugated
Full name:
interleukin 1 beta
Gene ID:
24494

Detection Antibody Information

Catalog Number:
84274-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241001H8
Reactivity:
rat
GenBank:
NM_031512

Conjugate:
Unconjugated
Full name:
interleukin 1 beta
Gene ID:
24494

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01210-3 targets IL-1 beta in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat IL-1 beta Recombinant antibody, PBS Only (Capture) 84274-2-PBS (241001G6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Rat IL-1 beta Recombinant antibody, PBS Only (Detector) 84274-1-PBS (241001H8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01210-3, RAT IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84274-2-PBS. Detection antibody: 84274-1-PBS. Standard: Eg31371. Range: 0.313-40 ng/mL.