

TSG101 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01221-1**

Capture Antibody Information

Catalog Number:
84334-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241657A11

Reactivity:
human

GenBank:
BC002487

Immunogen Catalog Number:
Ag28569

Conjugate:
Unconjugated

Full name:
tumor susceptibility gene 101

Gene ID:
7251

Detection Antibody Information

Catalog Number:
84334-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241657G12

Reactivity:
human

GenBank:
BC002487

Immunogen Catalog Number:
Ag28569

Conjugate:
Unconjugated

Full name:
tumor susceptibility gene 101

Gene ID:
7251

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01221-1 targets TSG101 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSG101 Recombinant antibody, PBS Only (Capture) 84334-2-PBS (241657A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: TSG101 Recombinant antibody, PBS Only (Detector) 84334-1-PBS (241657G12). 100 µg. Concentration 1 mg/mL.

Alternative TSG101 matched antibody pairs: MP01221-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01221-1, TSG101 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84334-2-PBS. Detection antibody: 84334-1-PBS. Standard: Ag28569. Range: 0.781-100 ng/mL.