

# MACROD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01269-1

## Capture Antibody Information

Catalog Number:  
84398-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
241506A3  
Reactivity:  
human  
GenBank:  
BC008316  
Immunogen Catalog Number:  
Ag34080

Conjugate:  
Unconjugated  
Full name:  
MACRO domain containing 1  
Gene ID:  
28992

## Detection Antibody Information

Catalog Number:  
84398-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
241506F12  
Reactivity:  
human  
GenBank:  
BC008316  
Immunogen Catalog Number:  
Ag34080

Conjugate:  
Unconjugated  
Full name:  
MACRO domain containing 1  
Gene ID:  
28992

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01269-1 targets MACROD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MACROD1 Recombinant antibody, PBS Only (Capture) 84398-1-PBS (241506A3). 100 µg. Concentration 1 mg/mL.

Detection antibody: MACROD1 Recombinant antibody, PBS Only (Detector) 84398-2-PBS (241506F12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01269-1, MACROD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84398-1-PBS. Detection antibody: 84398-2-PBS. Standard: Ag34080. Range: 0.781-100 ng/mL.