

TANC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01281-2**

Capture Antibody Information

Catalog Number:
84435-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241631A8
Reactivity:
human
GenBank:
NM_001394998
Immunogen Catalog Number:
Ag35427

Conjugate:
Unconjugated
Full name:
tetra-tryptophan repeat, ankyrin repeat and coiled-coil containing 2
Gene ID:
26115

Detection Antibody Information

Catalog Number:
84435-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241631E7
Reactivity:
human
GenBank:
NM_001394998
Immunogen Catalog Number:
Ag35427

Conjugate:
Unconjugated
Full name:
tetra-tryptophan repeat, ankyrin repeat and coiled-coil containing 2
Gene ID:
26115

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01281-2 targets TANC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TANC2 Recombinant antibody, PBS Only (Capture) 84435-1-PBS (241631A8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TANC2 Recombinant antibody, PBS Only (Detector) 84435-5-PBS (241631E7). 100 µg. Concentration 1 mg/mL.

Alternative TANC2 matched antibody pairs: MP01281-1, MP01281-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01281-2, TANC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84435-1-PBS. Detection antibody: 84435-5-PBS. Standard: Ag35427. Range: 0.781-100 ng/mL.