

Mouse TWEAKR/CD266 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01283-2**

Capture Antibody Information

Catalog Number:
84423-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241646C9
Reactivity:
mouse
GenBank:
NM_013749.2

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 12a
Gene ID:
27279

Detection Antibody Information

Catalog Number:
84423-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241646F11
Reactivity:
mouse
GenBank:
NM_013749.2

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 12a
Gene ID:
27279

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01283-2 targets TWEAKR/CD266 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse TWEAKR/CD266 Recombinant antibody, PBS Only (Capture) 84423-2-PBS (241646C9). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse TWEAKR/CD266 Recombinant antibody, PBS Only (Detector) 84423-5-PBS (241646F11). 100 µg. Concentration 1 mg/mL.

Alternative TWEAKR/CD266 matched antibody pairs: MP00861-1, MP01283-1, MP01283-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01283-2, MOUSE TWEAK R/TNFRSF12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84423-2-PBS. Detection antibody: 84423-5-PBS. Standard: Eg1838. Range: 0.156-20 ng/mL