

IL-1 alpha Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01289-2**

Capture Antibody Information

Catalog Number: 84411-3-PBS	Clone ID: 241752F11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 1, alpha
Isotype: IgG	GenBank: BC013142	Gene ID: 3552
Purification Method: Protein A purification	Immunogen Catalog Number: HZ-1320	

Detection Antibody Information

Catalog Number: 84411-4-PBS	Clone ID: 241752G3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 1, alpha
Isotype: IgG	GenBank: BC013142	Gene ID: 3552
Purification Method: Protein A purification	Immunogen Catalog Number: HZ-1320	

Applications

Tested Applications: Cytometric bead array	Range: 0.078-10 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01289-2 targets IL-1 alpha in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL1A Recombinant antibody, PBS Only (Capture/Detector) 84411-3-PBS (241752F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL1A Recombinant antibody, PBS Only (Detector) 84411-4-PBS (241752G3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01289-2, IL-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84411-3-PBS. Detection antibody: 84411-4-PBS. Standard: HZ-1320. Range: 0.078-10 ng/mL.