

Mouse Nectin-2/CD112

Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01292-3**

Capture Antibody Information

Catalog Number:
84439-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241798F12
Reactivity:
mouse
GenBank:
NM_008990.3

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 2
Gene ID:
19294

Detection Antibody Information

Catalog Number:
84439-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241798E6
Reactivity:
mouse
GenBank:
NM_008990.3

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 2
Gene ID:
19294

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01292-3 targets Nectin-2/CD112 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Capture) 84439-5-PBS (241798F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Detector) 84439-4-PBS (241798E6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01292-3, MOUSE Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84439-5-PBS. Detection antibody: 84439-4-PBS. Standard: Eg2019. Range: 0.313-20 ng/mL