

# CD40 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01294-1**

## Capture Antibody Information

**Catalog Number:**  
84406-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241760A4  
**Reactivity:**  
human  
**GenBank:**  
NM\_001250.6

**Conjugate:**  
Unconjugated  
**Full name:**  
CD40 molecule, TNF receptor superfamily member 5  
**Gene ID:**  
958

## Detection Antibody Information

**Catalog Number:**  
84406-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241760A11  
**Reactivity:**  
human  
**GenBank:**  
NM\_001250.6

**Conjugate:**  
Unconjugated  
**Full name:**  
CD40 molecule, TNF receptor superfamily member 5  
**Gene ID:**  
958

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01294-1 targets CD40 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD40 Recombinant antibody, PBS Only (Capture) 84406-1-PBS (241760A4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD40 Recombinant antibody, PBS Only (Detector) 84406-2-PBS (241760A11). 100 µg. Concentration 1 mg/mL.

Alternative CD40 matched antibody pairs: MP01294-2, MP01294-3, MP50064-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01294-1, CD40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84406-1-PBS. Detection antibody: 84406-2-PBS. Standard: Eg1797. Range: 0.625-80 ng/mL