

# Mouse CEACAM1/CD66a Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01297-3**

## Capture Antibody Information

**Catalog Number:**  
84409-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241856E3  
**Reactivity:**  
mouse  
**GenBank:**  
**Immunogen Catalog Number:**  
Eg1423

**Conjugate:**  
Unconjugated  
**Full name:**  
carcinoembryonic antigen-related  
cell adhesion molecule 1  
**Gene ID:**  
26365

## Detection Antibody Information

**Catalog Number:**  
84409-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241856H5  
**Reactivity:**  
mouse  
**GenBank:**  
**Immunogen Catalog Number:**  
Eg1423

**Conjugate:**  
Unconjugated  
**Full name:**  
carcinoembryonic antigen-related  
cell adhesion molecule 1  
**Gene ID:**  
26365

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead  
Array)

**Recommended Dilutions:**  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP01297-3 targets CEACAM1/CD66a in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CEACAM-1/CD66a Recombinant antibody, PBS Only 84409-2-PBS (241856E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CEACAM1/CD66a Recombinant antibody, PBS Only (Capture) 84409-4-PBS (241856H5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01297-3, MOUSE CEACAM1/CD66a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84409-2-PBS. Detection antibody: 84409-4-PBS. Standard: Eg1423. Range: 0.313-40 ng/mL