

Mouse CD157 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01303-2**

Capture Antibody Information

Catalog Number:
84417-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241822A6
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
bone marrow stromal cell antigen 1
Gene ID:
12182

Detection Antibody Information

Catalog Number:
84417-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241822H1
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
bone marrow stromal cell antigen 1
Gene ID:
12182

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01303-2 targets CD157 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD157 Recombinant antibody, PBS Only (Capture) 84417-2-PBS (241822A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD157 Recombinant antibody, PBS Only (Capture/Detector) 84417-4-PBS (241822H1). 100 µg. Concentration 1 mg/mL.

Alternative CD157 matched antibody pairs: MP01225-1, MP01225-2, MP01225-3, MP01225-4, MP01303-1, MP01303-3, MP01303-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01303-2, MOUSE CD157 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84417-2-PBS. Detection antibody: 84417-4-PBS. Standard: Eg1528. Range: 0.313-40 ng/mL.