

CKAP4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01324-1**

Capture Antibody Information

Catalog Number: 84485-1-PBS	Clone ID: 241871A11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cytoskeleton-associated protein 4
Isotype: IgG	GenBank: BC082972	Gene ID: 10970
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10022	

Detection Antibody Information

Catalog Number: 84485-2-PBS	Clone ID: 241871F12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cytoskeleton-associated protein 4
Isotype: IgG	GenBank: BC082972	Gene ID: 10970
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10022	

Applications

Tested Applications: Cytometric bead array, Sandwich ELISA	Range: 0.156-20 ng/mL (Cytometric Bead Array) 9.77-313 pg/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01324-1 targets CKAP4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: CKAP4 Recombinant antibody, PBS Only (Capture) 84485-1-PBS (241871A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CKAP4 Recombinant antibody, PBS Only (Detector) 84485-2-PBS (241871F12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only


Selected Validation Data



Sandwich ELISA standard curve of MP01324-1, Human CKAP4 Recombinant Matched Antibody Pair - PBS only. 84485-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10022. 84485-2-PBS was HRP conjugated as the detection antibody. Range: 9.77-313 pg/mL.



The mean CKAP4 concentration was determined to be 13.6 ng/mL in A431 cell extract based on a 1.3 mg/mL extract load and 9.2 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP01324-1, CKAP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84485-1-PBS. Detection antibody: 84485-2-PBS. Standard: Ag10022. Range: 0.156-20 ng/mL.