

For Research Use Only

Mouse NKp46/NCR1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number: **MP01343-3**

Capture Antibody Information

Catalog Number:
84479-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241804A12
Reactivity:
mouse
GenBank:
NM_010746.3

Conjugate:
Unconjugated
Full name:
natural cytotoxicity triggering receptor 1
Gene ID:
17086

Detection Antibody Information

Catalog Number:
84479-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241804F1
Reactivity:
mouse
GenBank:
NM_010746.3

Conjugate:
Unconjugated
Full name:
natural cytotoxicity triggering receptor 1
Gene ID:
17086

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01343-3 targets NKp46/NCR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse NKp46/NCR1 Recombinant antibody, PBS Only (Capture) 84479-1-PBS (241804A12). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse NKp46/NCR1 Recombinant antibody, PBS Only (Detector) 84479-3-PBS (241804F1). 100 µg. Concentration 1 mg/mL.

Alternative NKp46/NCR1 matched antibody pairs: MP01299-1, MP01299-2, MP01299-3, MP01343-1, MP01343-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01343-3, MOUSE Nkp46/NCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84479-1-PBS. Detection antibody: 84479-3-PBS. Standard: Eg1825. Range: 0.625-40 ng/mL