

# Mouse Neutrophil elastase Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01355-3**

## Capture Antibody Information

**Catalog Number:**  
84473-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241842E7  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_015779.2

**Conjugate:**  
Unconjugated  
**Full name:**  
elastase 2, neutrophil  
**Gene ID:**  
50701

## Detection Antibody Information

**Catalog Number:**  
84473-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241842E3  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_015779.2

**Conjugate:**  
Unconjugated  
**Full name:**  
elastase 2, neutrophil  
**Gene ID:**  
50701

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01355-3 targets Neutrophil elastase in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Neutrophil elastase Recombinant antibody, PBS Only (Capture/Detector) 84473-4-PBS (241842E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Neutrophil elastase Recombinant antibody, PBS Only (Detector) 84473-3-PBS (241842E3). 100 µg. Concentration 1 mg/mL.

Alternative Neutrophil elastase matched antibody pairs: MP01355-1, MP01355-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01355-3, MOUSE Neutrophil elastase Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84473-4-PBS. Detection antibody: 84473-3-PBS. Standard: Eg2126. Range: 0.781-100 ng/mL