

TRAF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01404-1**

Capture Antibody Information

Catalog Number: 84581-2-PBS	Clone ID: 241899G11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: TNF receptor-associated factor 2
Isotype: IgG	GenBank: BC032410	Gene ID: 7186
Purification Method: Protein A purification	Immunogen Catalog Number: Ag25387	

Detection Antibody Information

Catalog Number: 84581-1-PBS	Clone ID: 241899E6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: TNF receptor-associated factor 2
Isotype: IgG	GenBank: BC032410	Gene ID: 7186
Purification Method: Protein A purification	Immunogen Catalog Number: Ag25387	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01404-1 targets TRAF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRAF2 Recombinant antibody, PBS Only 84581-2-PBS (241899G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRAF2 Recombinant antibody, PBS Only (Detector) 84581-1-PBS (241899E6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01404-1, TRAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84581-2-PBS. Detection antibody: 84581-1-PBS. Standard: Ag25387. Range: 1.563-200 ng/mL.