

Mouse TLR3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01424-2**

Capture Antibody Information

Catalog Number:
84618-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241792A5
Reactivity:
mouse
GenBank:
NM_126166.4

Conjugate:
Unconjugated
Full name:
toll-like receptor 3
Gene ID:
142980

Detection Antibody Information

Catalog Number:
84618-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241792F9
Reactivity:
mouse
GenBank:
NM_126166.4

Conjugate:
Unconjugated
Full name:
toll-like receptor 3
Gene ID:
142980

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01424-2 targets TLR3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse TLR3 Recombinant antibody, PBS Only (Capture) 84618-1-PBS (241792A5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse TLR3 Recombinant antibody, PBS Only (Detector) 84618-3-PBS (241792F9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01424-2, MOUSE TLR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84618-1-PBS. Detection antibody: 84618-3-PBS. Standard: Eg1803. Range: 0.156-20 ng/mL.