

# EDDM3B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01431-1**

## Capture Antibody Information

**Catalog Number:**  
84597-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241998B7

**Reactivity:**  
human

**GenBank:**  
NM\_022360.5

**Conjugate:**  
Unconjugated

**Full name:**  
family with sequence similarity 12, member B (epididymal)

**Gene ID:**  
64184

## Detection Antibody Information

**Catalog Number:**  
84597-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241998B5

**Reactivity:**  
human

**GenBank:**  
NM\_022360.5

**Conjugate:**  
Unconjugated

**Full name:**  
family with sequence similarity 12, member B (epididymal)

**Gene ID:**  
64184

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.469-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01431-1 targets EDDM3B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EDDM3B Recombinant antibody, PBS Only (Capture) 84597-2-PBS (241998B7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: EDDM3B Recombinant antibody, PBS Only (Detector) 84597-1-PBS (241998B5). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative EDDM3B matched antibody pairs: MP01431-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01431-1, EDDM3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84597-2-PBS. Detection antibody: 84597-1-PBS. Standard: Eg2443. Range: 0.469-60 ng/mL