

KAT2A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01463-1**

Capture Antibody Information

Catalog Number: 84669-2-PBS	Clone ID: 242113C7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, rat	Full name: K(lysine) acetyltransferase 2A
Isotype: IgG	GenBank: BC032743	Gene ID: 2648
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29005	

Detection Antibody Information

Catalog Number: 84669-3-PBS	Clone ID: 242113H2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: K(lysine) acetyltransferase 2A
Isotype: IgG	GenBank: BC032743	Gene ID: 2648
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29005	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01463-1 targets KAT2A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KAT2A/GCN5 Recombinant antibody, PBS Only (Capture) 84669-2-PBS (242113C7). 100 µg. Concentration 1 mg/mL.

Detection antibody: KAT2A Recombinant antibody, PBS Only (Detector) 84669-3-PBS (242113H2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01463-1, KAT2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84669-2-PBS. Detection antibody: 84669-3-PBS. Standard:Ag29005. Range: 0.313-40 ng/mL