

EPCR/CD201 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01480-2**

Capture Antibody Information

Catalog Number:

84648-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241806D1

Reactivity:

human

GenBank:

BC014451

Conjugate:

Unconjugated

Full name:

protein C receptor, endothelial (EPCR)

Gene ID:

10544

Detection Antibody Information

Catalog Number:

84648-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241806E10

Reactivity:

human

GenBank:

BC014451

Conjugate:

Unconjugated

Full name:

protein C receptor, endothelial (EPCR)

Gene ID:

10544

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01480-2 targets EPCR/CD201 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPCR/CD201 Recombinant antibody, PBS Only (Capture) 84648-2-PBS (241806D1). 100 µg. Concentration 1 mg/mL.

Detection antibody: EPCR/CD201 Recombinant antibody, PBS Only (Detector) 84648-4-PBS (241806E10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01480-2, PROCR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84648-2-PBS. Detection antibody: 84648-4-PBS. Standard: Eg2339. Range: 0.625-80 ng/mL