

CDC40 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01501-2**

Capture Antibody Information

Catalog Number:
84679-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242034H8
Reactivity:
human
GenBank:
BC126114
Immunogen Catalog Number:
Ag36151

Conjugate:
Unconjugated
Full name:
cell division cycle 40 homolog (S. cerevisiae)
Gene ID:
51362

Detection Antibody Information

Catalog Number:
84679-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242034G6
Reactivity:
human
GenBank:
BC126114
Immunogen Catalog Number:
Ag36151

Conjugate:
Unconjugated
Full name:
cell division cycle 40 homolog (S. cerevisiae)
Gene ID:
51362

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01501-2 targets CDC40 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDC40 Recombinant antibody, PBS Only (Capture) 84679-3-PBS (242034H8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CDC40 Recombinant antibody, PBS Only (Detector) 84679-1-PBS (242034G6). 100 µg. Concentration 1 mg/mL.

Alternative CDC40 matched antibody pairs: MP01501-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01501-2, CDC40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84679-3-PBS. Detection antibody: 84679-1-PBS. Standard: Ag36151. Range: 0.313-40 ng/mL