

# Mouse MECR Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01508-1**

## Capture Antibody Information

**Catalog Number:**  
84729-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242151D11  
**Reactivity:**  
mouse  
**GenBank:**  
BC003864  
**Immunogen Catalog Number:**  
Ag0465

**Conjugate:**  
Unconjugated  
**Full name:**  
mitochondrial trans-2-enoyl-CoA reductase  
**Gene ID:**  
26922

## Detection Antibody Information

**Catalog Number:**  
84729-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242151G10  
**Reactivity:**  
mouse  
**GenBank:**  
BC003864  
**Immunogen Catalog Number:**  
Ag0465

**Conjugate:**  
Unconjugated  
**Full name:**  
mitochondrial trans-2-enoyl-CoA reductase  
**Gene ID:**  
26922

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01508-1 targets MECR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse MECR Recombinant antibody, PBS Only (Capture) 84729-2-PBS (242151D11). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: Mouse MECR Recombinant antibody, PBS Only (Detector) 84729-3-PBS (242151G10). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01508-1, MOUSE MECR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84729-2-PBS. Detection antibody: 84729-3-PBS. Standard: Ag0465. Range: 0.625-80 ng/mL