

# C11orf68 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01522-1

## Capture Antibody Information

**Catalog Number:**  
84725-3-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
242221G9

**Reactivity:**  
human

**GenBank:**  
BC010512

**Immunogen Catalog Number:**  
Ag36023

**Conjugate:**  
Unconjugated

**Full name:**  
chromosome 11 open reading frame 68

**Gene ID:**  
83638

## Detection Antibody Information

**Catalog Number:**  
84725-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
242221A1

**Reactivity:**  
human

**GenBank:**  
BC010512

**Immunogen Catalog Number:**  
Ag36023

**Conjugate:**  
Unconjugated

**Full name:**  
chromosome 11 open reading frame 68

**Gene ID:**  
83638

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01522-1 targets C11orf68 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C11orf68 Recombinant antibody, PBS Only (Capture) 84725-3-PBS (242221G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: C11orf68 Recombinant antibody, PBS Only (Detector) 84725-1-PBS (242221A1). 100 µg. Concentration 1 mg/mL.

Alternative C11orf68 matched antibody pairs: MP01522-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01522-1, C11orf68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84725-3-PBS. Detection antibody: 84725-1-PBS. Standard: Ag36023. Range: 0.781-100 ng/mL