

# C11orf68 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01522-2**

## Capture Antibody Information

Catalog Number:  
84725-3-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
242221G9  
Reactivity:  
human  
GenBank:  
BC010512  
Immunogen Catalog Number:  
Ag36023

Conjugate:  
Unconjugated  
Full name:  
chromosome 11 open reading frame 68  
Gene ID:  
83638

## Detection Antibody Information

Catalog Number:  
84725-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
242221A4  
Reactivity:  
human  
GenBank:  
BC010512  
Immunogen Catalog Number:  
Ag36023

Conjugate:  
Unconjugated  
Full name:  
chromosome 11 open reading frame 68  
Gene ID:  
83638

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01522-2 targets C11orf68 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C11orf68 Recombinant antibody, PBS Only (Capture) 84725-3-PBS (242221G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: C11orf68 Recombinant antibody, PBS Only (Detector) 84725-2-PBS (242221A4). 100 µg. Concentration 1 mg/mL.

Alternative C11orf68 matched antibody pairs: MP01522-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01522-2, C11orf68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84725-3-PBS. Detection antibody: 84725-2-PBS. Standard: Ag36023. Range: 0.781-100 ng/mL.