

PRDM8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01532-2**

Capture Antibody Information

Catalog Number: 84716-1-PBS	Clone ID: 242203B4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: PR domain containing 8
Isotype: IgG	GenBank: BC071584	Gene ID: 56978
Purification Method: Protein A purification	Immunogen Catalog Number: Ag23373	

Detection Antibody Information

Catalog Number: 84716-3-PBS	Clone ID: 242203F2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: PR domain containing 8
Isotype: IgG	GenBank: BC071584	Gene ID: 56978
Purification Method: Protein A purification	Immunogen Catalog Number: Ag23373	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP01532-2 targets PRDM8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRDM8 Recombinant antibody, PBS Only (Capture) 84716-1-PBS (242203B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: PRDM8 Recombinant antibody, PBS Only (Detector) 84716-3-PBS (242203F2). 100 µg. Concentration 1 mg/mL.

Alternative PRDM8 matched antibody pairs: MP01532-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01532-2, PRDM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84716-1-PBS. Detection antibody: 84716-3-PBS. Standard: Ag23373. Range: 0.781-100 ng/mL.