

# DR5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01541-2**

## Capture Antibody Information

**Catalog Number:**  
84784-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
242319G5

**Reactivity:**  
human

**GenBank:**  
BC001281

**Immunogen Catalog Number:**  
Eg1648

**Conjugate:**  
Unconjugated

**Full name:**  
tumor necrosis factor receptor superfamily, member 10b

**Gene ID:**  
8795

## Detection Antibody Information

**Catalog Number:**  
84784-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
242319B12

**Reactivity:**  
human

**GenBank:**  
BC001281

**Immunogen Catalog Number:**  
Eg1648

**Conjugate:**  
Unconjugated

**Full name:**  
tumor necrosis factor receptor superfamily, member 10b

**Gene ID:**  
8795

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01541-2 targets DR5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DR5 Recombinant antibody, PBS Only 84784-2-PBS (242319G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: DR5 Recombinant antibody, PBS Only (Detector) 84784-1-PBS (242319B12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01541-2, DR5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84784-2-PBS. Detection antibody: 84784-1-PBS. Standard: Eg1648. Range: 0.625-40 ng/mL