

SPTAN1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01560-2**

Capture Antibody Information

Catalog Number:
84761-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242104B7

Reactivity:
human, mouse, rat

GenBank:
NM_003127.3

Immunogen Catalog Number:
Ag36472

Conjugate:
Unconjugated

Full name:
spectrin, alpha, non-erythrocytic 1 (alpha-fodrin)

Gene ID:
6709

Detection Antibody Information

Catalog Number:
84761-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242104D6

Reactivity:
human

GenBank:
NM_003127.3

Immunogen Catalog Number:
Ag36472

Conjugate:
Unconjugated

Full name:
spectrin, alpha, non-erythrocytic 1 (alpha-fodrin)

Gene ID:
6709

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01560-2 targets SPTAN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPTAN1 Recombinant antibody, PBS Only 84761-2-PBS (242104B7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SPTAN1 Recombinant antibody, PBS Only (Detector) 84761-3-PBS (242104D6). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01560-2, SPTAN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84761-2-PBS. Detection antibody: 84761-3-PBS. Standard: Ag36472. Range: 0.469-60 ng/mL.