

GPD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01579-2**

Capture Antibody Information

Catalog Number:
84800-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242359E12

Reactivity:
human

GenBank:
BC032234

Immunogen Catalog Number:
Ag4278

Conjugate:
Unconjugated

Full name:
glycerol-3-phosphate dehydrogenase 1 (soluble)

Gene ID:
2819

Detection Antibody Information

Catalog Number:
84800-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242359E2

Reactivity:
human

GenBank:
BC032234

Immunogen Catalog Number:
Ag4278

Conjugate:
Unconjugated

Full name:
glycerol-3-phosphate dehydrogenase 1 (soluble)

Gene ID:
2819

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01579-2 targets GPD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GPD1 Recombinant antibody, PBS Only (Capture/Detector) 84800-2-PBS (242359E12). 100 µg. Concentration 1 mg/mL.

Detection antibody: GPD1 Recombinant antibody, PBS Only (Detector) 84800-3-PBS (242359E2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01579-2, GPD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84800-2-PBS. Detection antibody: 84800-3-PBS. Standard: Ag4278. Range: 1.563-100 ng/mL.