

RNF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01591-1**

Capture Antibody Information

Catalog Number: 84793-3-PBS	Clone ID: 242181D2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ring finger protein 2
Isotype: IgG	GenBank: BC012583	Gene ID: 6045
Purification Method: Protein A purification	Immunogen Catalog Number: Ag8900	

Detection Antibody Information

Catalog Number: 84793-1-PBS	Clone ID: 242181B7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ring finger protein 2
Isotype: IgG	GenBank: BC012583	Gene ID: 6045
Purification Method: Protein A purification	Immunogen Catalog Number: Ag8900	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01591-1 targets RNF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RNF2 Recombinant antibody, PBS Only 84793-3-PBS (242181D2). 100 µg. Concentration 1 mg/mL.

Detection antibody: RNF2 Recombinant antibody, PBS Only 84793-1-PBS (242181B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01591-1, RNF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84793-3-PBS. Detection antibody: 84793-1-PBS. Standard: Ag8900. Range: 0.313-40 ng/mL