

SHOC2/Sur-8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01674-2

Capture Antibody Information

Catalog Number: 84917-2-PBS	Clone ID: 242343D10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: soc-2 suppressor of clear homolog (C. elegans)
Isotype: IgG	GenBank: BC050445	Gene ID: 8036
Purification Method: Protein A purification	Immunogen Catalog Number: Ag11754	

Detection Antibody Information

Catalog Number: 84917-4-PBS	Clone ID: 242343F2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: soc-2 suppressor of clear homolog (C. elegans)
Isotype: IgG	GenBank: BC050445	Gene ID: 8036
Purification Method: Protein A purification	Immunogen Catalog Number: Ag11754	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01674-2 targets SHOC2/Sur-8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SHOC2/Sur-8 Recombinant antibody, PBS Only (Capture) 84917-2-PBS (242343D10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SHOC2 Recombinant antibody, PBS Only (Detector) 84917-4-PBS (242343F2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01674-2, SHOC2/Sur-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84917-2-PBS. Detection antibody: 84917-4-PBS. Standard: Ag11754. Range: 0.781-100 ng/mL.