

For Research Use Only

Mouse CSF1R/CD115 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01676-1

Capture Antibody Information

Catalog Number:

84919-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242213F5

Reactivity:

mouse

GenBank:

NM_001037859.2

Conjugate:

Unconjugated

Full name:

colony stimulating factor 1 receptor

Gene ID:

12978

Detection Antibody Information

Catalog Number:

84919-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242213D2

Reactivity:

mouse

GenBank:

NM_001037859.2

Conjugate:

Unconjugated

Full name:

colony stimulating factor 1 receptor

Gene ID:

12978

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01676-1 targets CSF1R/CD115 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse M-CSF R/CD115 Recombinant antibody, PBS Only (Capture) 84919-3-PBS (242213F5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Mouse M-CSF R/CD115 Recombinant antibody, PBS Only (Detector) 84919-2-PBS (242213D2). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01676-1, MOUSE CSF1R/CD115 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84919-3-PBS. Detection antibody: 84919-2-PBS. Standard: Eg1816. Range: 0.625-80 ng/mL