

# SEC31A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01696-2**

## Capture Antibody Information

Catalog Number: 84942-3-PBS	Clone ID: 242531G6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SEC31 homolog A (S. cerevisiae)
Isotype: IgG	GenBank: BC084583	Gene ID: 22872
Purification Method: Protein A purification	Immunogen Catalog Number: Ag12314	

## Detection Antibody Information

Catalog Number: 84942-1-PBS	Clone ID: 242531A11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SEC31 homolog A (S. cerevisiae)
Isotype: IgG	GenBank: BC084583	Gene ID: 22872
Purification Method: Protein A purification	Immunogen Catalog Number: Ag12314	

## Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

## Product Information

MP01696-2 targets SEC31A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEC31A Recombinant antibody, PBS Only (Capture) 84942-3-PBS (242531G6). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: SEC31A Recombinant antibody, PBS Only (Detector) 84942-1-PBS (242531A11). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01696-2, SEC31A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84942-3-PBS. Detection antibody: 84942-1-PBS. Standard: Ag12314. Range: 0.313-40 ng/mL.