

Selenoprotein P Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01743-2**

Capture Antibody Information

Catalog Number: 85017-3-PBS	Clone ID: 242510H7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: selenoprotein P, plasma, 1
Isotype: IgG	GenBank: NM_001085486.2	Gene ID: 6414
Purification Method: Protein A purification	Immunogen Catalog Number: Ag32952	

Detection Antibody Information

Catalog Number: 85017-2-PBS	Clone ID: 242510D6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: selenoprotein P, plasma, 1
Isotype: IgG	GenBank: NM_001085486.2	Gene ID: 6414
Purification Method: Protein A purification	Immunogen Catalog Number: Ag32952	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01743-2 targets Selenoprotein P in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Selenoprotein P Recombinant antibody, PBS Only (Capture) 85017-3-PBS (242510H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Selenoprotein P Recombinant antibody, PBS Only (Capture/Detector) 85017-2-PBS (242510D6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01743-2, Selenoprotein P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85017-3-PBS. Detection antibody: 85017-2-PBS. Standard: Ag32952. Range: 0.781-100 ng/mL