

SMARCA5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01767-2**

Capture Antibody Information

Catalog Number:
85019-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242512A8

Reactivity:
human

GenBank:
BC023144

Immunogen Catalog Number:
Ag3759

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 5

Gene ID:
8467

Detection Antibody Information

Catalog Number:
85019-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242512A10

Reactivity:
human

GenBank:
BC023144

Immunogen Catalog Number:
Ag3759

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 5

Gene ID:
8467

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01767-2 targets SMARCA5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA5 Recombinant antibody, PBS Only (Capture) 85019-3-PBS (242512A8). 100 µg. Concentration 1 mg/mL.

Detection antibody: SMARCA5 Recombinant antibody, PBS Only (Detector) 85019-1-PBS (242512A10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01767-2, SMARCA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85019-3-PBS. Detection antibody: 85019-1-PBS. Standard: Ag3759. Range: 0.781-100 ng/mL.