

SEC22B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01774-1**

Capture Antibody Information

Catalog Number:
85037-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242659E4
Reactivity:
human
GenBank:
BC001364
Immunogen Catalog Number:
Ag12981

Conjugate:
Unconjugated
Full name:
SEC22 vesicle trafficking protein homolog B (S. cerevisiae)
Gene ID:
9554

Detection Antibody Information

Catalog Number:
85037-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242659F8
Reactivity:
human
GenBank:
BC001364
Immunogen Catalog Number:
Ag12981

Conjugate:
Unconjugated
Full name:
SEC22 vesicle trafficking protein homolog B (S. cerevisiae)
Gene ID:
9554

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01774-1 targets SEC22B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEC22B Recombinant antibody, PBS Only (Capture) 85037-3-PBS (242659E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: SEC22B Recombinant antibody, PBS Only (Detector) 85037-4-PBS (242659F8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01774-1, SEC22B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85037-3-PBS. Detection antibody: 85037-4-PBS. Standard: Ag12981. Range: 0.313-40 ng/mL