

Mouse Fetuin A/AHSG Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01802-2**

Capture Antibody Information

Catalog Number:
85085-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242732A11

Reactivity:
mouse

GenBank:
NM_013465.2

Conjugate:
Unconjugated

Full name:
alpha-2-HS-glycoprotein

Gene ID:
11625

Detection Antibody Information

Catalog Number:
85085-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242732H3

Reactivity:
mouse

GenBank:
NM_013465.2

Conjugate:
Unconjugated

Full name:
alpha-2-HS-glycoprotein

Gene ID:
11625

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01802-2 targets Fetuin A/AHSG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Fetuin A/AHSG Recombinant antibody, PBS Only (Capture) 85085-1-PBS (242732A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Fetuin A/AHSG Recombinant antibody, PBS Only (Detector) 85085-3-PBS (242732H3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01802-2, MOUSE Ahsg Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85085-1-PBS. Detection antibody: 85085-3-PBS. Standard: Eg3151. Range: 0.313-40 ng/mL