

# CCNL2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01876-1

## Capture Antibody Information

Catalog Number: 85149-3-PBS	Clone ID: 242806D8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cyclin L2
Isotype: IgG	GenBank: BC016333	Gene ID: 81669
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1710	

## Detection Antibody Information

Catalog Number: 85149-1-PBS	Clone ID: 242806B12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cyclin L2
Isotype: IgG	GenBank: BC016333	Gene ID: 81669
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1710	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP01876-1 targets CCNL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCNL2 Recombinant antibody, PBS Only (Capture) 85149-3-PBS (242806D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCNL2 Recombinant antibody, PBS Only (Detector) 85149-1-PBS (242806B12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01876-1, CCNL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85149-3-PBS. Detection antibody: 85149-1-PBS. Standard: Ag1710. Range: 0.781-100 ng/mL.