

HSCB Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01886-2**

Capture Antibody Information

Catalog Number: 85185-2-PBS	Clone ID: 242836A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HscB iron-sulfur cluster co-chaperone homolog (E. coli)
Isotype: IgG	GenBank: BC065569	Gene ID: 150274
Purification Method: Protein A purification	Immunogen Catalog Number: Ag7268	

Detection Antibody Information

Catalog Number: 85185-1-PBS	Clone ID: 242836A5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HscB iron-sulfur cluster co-chaperone homolog (E. coli)
Isotype: IgG	GenBank: BC065569	Gene ID: 150274
Purification Method: Protein A purification	Immunogen Catalog Number: Ag7268	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP01886-2 targets HSCB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSCB Recombinant antibody, PBS Only (Capture) 85185-2-PBS (242836A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: HSCB Recombinant antibody, PBS Only (Detector) 85185-1-PBS (242836A5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01886-2, HSCB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85185-2-PBS. Detection antibody: 85185-1-PBS. Standard: Ag7268. Range: 0.313-40 ng/mL