

Mouse CLEC12A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01887-2**

Capture Antibody Information

Catalog Number:
85176-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242620C12
Reactivity:
mouse
GenBank:
/

Conjugate:
Unconjugated
Full name:
C-type lectin domain family 12,
member a
Gene ID:
232413

Detection Antibody Information

Catalog Number:
85176-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242620D1
Reactivity:
mouse
GenBank:
/

Conjugate:
Unconjugated
Full name:
C-type lectin domain family 12,
member a
Gene ID:
232413

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01887-2 targets CLEC12A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Clec12a Recombinant antibody, PBS Only (Capture) 85176-1-PBS (242620C12). 100 µg.
Concentration 1 mg/mL.

Detection antibody: Mouse Clec12a Recombinant antibody, PBS Only (Detector) 85176-2-PBS (242620D1). 100 µg.
Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01887-2, MOUSE CLEC12A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85176-1-PBS. Detection antibody: 85176-2-PBS. Standard: Eg2726. Range: 0.313-40 ng/mL