

For Research Use Only

Mouse CD72 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01912-1

Capture Antibody Information

Catalog Number:
85208-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242846E3

Reactivity:
mouse

GenBank:
NM_007654.2

Immunogen Catalog Number:
Eg2708

Conjugate:
Unconjugated

Full name:
CD72 antigen

Gene ID:
12517

Detection Antibody Information

Catalog Number:
85208-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242846D10

Reactivity:
mouse

GenBank:
NM_007654.2

Immunogen Catalog Number:
Eg2708

Conjugate:
Unconjugated

Full name:
CD72 antigen

Gene ID:
12517

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01912-1 targets CD72 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Cd72 Recombinant monoclonal antibody, PBS Only (Capture) 85208-3-PBS (242846E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Cd72 Recombinant monoclonal antibody, PBS Only (Detector) 85208-2-PBS (242846D10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01912-1, MOUSE CD72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85208-3-PBS. Detection antibody: 85208-2-PBS. Standard: Eg2708. Range: 0.781-100 ng/mL.

The mean CD72 concentration was determined to be 11.7 ng/mL in RAW264.7 cell extract based on a 1.5 mg/mL extract load. Mouse brain tissue extract were served as a negative control.