

CHT1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01923-2**

Capture Antibody Information

Catalog Number: 85405-2-PBS	Clone ID: 242904G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: solute carrier family 5 (choline transporter), member 7
Isotype: IgG	GenBank: BC111525	Gene ID: 60482
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16136	

Detection Antibody Information

Catalog Number: 85405-3-PBS	Clone ID: 242904G5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: solute carrier family 5 (choline transporter), member 7
Isotype: IgG	GenBank: BC111525	Gene ID: 60482
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16136	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01923-2 targets CHT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC5A7 Recombinant antibody, PBS Only (Capture) 85405-2-PBS (242904G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC5A7 Recombinant antibody, PBS Only (Detector) 85405-3-PBS (242904G5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01923-2, CHT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85405-2-PBS. Detection antibody: 85405-3-PBS. Standard: Ag16136. Range: 0.781-100 ng/mL.